

## 1M in tetrahydrofuran

#### **OTHER NAMES**

#### **TYPICAL PROPERTIES**

Tetrahydrofuran borane

CAS REG. NUMBER


[14044-65-6]

### **EINECS**

237-881-8 for BTHF 203-726-8 for THF

## FORMULA

 $(CH_2)_4OBH_3$ 



## FORMULA WEIGHT

85.94 g/mol

### **RELATED REAGENTS**

Dimethylsulfide borane *N*,*N*-Diethylaniline borane

## PHYSICAL DATA Density at 25 °C:

Flash point (closed cup): Boiling point: 0.876 g/mL -5.7 °F (-22 °C) 66 °C for THF, see Stability sec.

### STABILITY

BTHF is thermally unstable and must be kept cold. BTHF stabilized with <0.005 M NaBH<sub>4</sub> per mole of "BH<sub>3</sub>". BTHF is susceptible to hydrolysis, readily reacting with water to form hydrogen and boric acid. It reacts readily with atmospheric moisture upon exposure to air resulting in a decrease in purity.

In the absence of a substrate, BTHF decomposes by cleavage of the ether ring and above 50  $^{\circ}\text{C}$  can evolve diborane.

## **APPEARANCE**

Colorless liquid

## **STORAGE CONDITIONS**

To maximize shelf-life, BTHF solutions should be stored below 5  $^{\circ}\text{C}.$ 

## PACKAGING

Packaged in cylinders

- 800 ml, containing 0.7 Kg 1M BTHF solution
- 18 liter, containing 16 Kgs 1M BTHF solution
- 90 liter, containing 79 Kgs 1M BTHF solution
- 400 liter, containing 351 Kgs 1M BTHF solution

### SHIPPING INFORMATION

UN-3148, PG I

1420 Mars-Evans City Road Evans City, PA 16033 USA Tel: 1-866-4BORANE or 412-967-4141 Fax: 412-967-4140 Email: <u>Callery.info@MSAnet.com</u> http://www.Callery.com

## PRODUCT FEATURES

- Reduces selected functional groups
- Hydroborating agent
- Borane source for oxazaborolidine catalyzed asymmetric reductions
- High purity
- Consistent quality
- Non-corrosive to metals
- Custom packaging available

### **PRODUCT BENEFITS**


- Excellent reagent for the reduction of amides to amines
- Preferred reagent for reduction of carboxylic acids to alcohols

BORANE COMPLEX

#### **APPLICATIONS**

Borane-tetrahydrofuran complex (BTHF) is a valuable reagent for the reduction of functional groups and for hydroboration reactions with carbon-carbon double and triple bonds. Functional groups that are readily reduced by BTHF include aldehyde, ketone, carboxylic acid, amide, oxime, imine, and nitrile. The carboxylic acid group is reduced at a faster rate than most groups including non-conjugated alkene.<sup>1</sup> Conjugated  $\alpha,\beta$ -unsaturated carboxylic acids give saturated alcohols as the major products.

Ketones and the carbonyl of enones are effectively reduced with borane-tetrahydrofuran. The addition of borohydride to the reaction solution is advantageous for accelerated reduction<sup>2</sup> as well as higher selectivity towards carbonyl reduction in conjugated and non-conjugated enones.<sup>3</sup>




Asymmetric ketone reduction using chiral oxazaborolidine catalysts was recently reviewed.<sup>4</sup> Work at Callery with BTHF improved on reaction conditions to provide consistent results in the reduction.<sup>5</sup>

# (*R*)-MeCBS reduction of acetophenone with NaBH<sub>4</sub> stabilized 1M BTHF under optimized conditions

entry	conditions	% ee°
1	Acetophenone and BTHF same rate <sup>a</sup>	95.2
2	Acetophenone 2X faster than BTHF	96.3
3	BTHF 2X faster than acetophenone	95.6
4	BF <sub>3</sub> Additive <sup>b</sup>	92.8


<sup>a.</sup> Acetophenone and stabilized BTHF (~0.005M NaBH<sub>4</sub>) added simultaneously to (*R*)-MeCBS. <sup>b.</sup> Stabilized BTHF solution was added to a mixture of ketone, BF<sub>3</sub> and (*R*)-MeCBS. <sup>c</sup> Enantiomeric excess was measured with chiral column J & W 30m x 0.25 mm CDX-B.

EPIX Medical<sup>6</sup> has used BTHF for the reduction of an amide in the synthesis of an injectable vascular contrast agent. They optimized conditions to minimize reduction of the BOC protecting group.



Several neutral work-up methods<sup>7</sup> are available in addition to acidic quenching of the amine borane intermediate produced.

Alkylboranes obtained via hydroboration are thermally sensitive and are prone to undergo dehydroboration-rehydroboration until the boron group is at the terminal position.<sup>8</sup> This migration has been used as an approach for regio- and stereoselective control of up to three chiral centers using a new stereoselective migration of organoboranes.<sup>9</sup>



<sup>1</sup>Brown, H.C. *et al. J. Org. Chem.* **1973**, 38, 2786.

<sup>2</sup> Jockel, H.; Schmidt, R. J. Chem. Soc. Perkin Trans. 2 1997, 2719.

- <sup>3</sup> Arase, A.; Hoshi, M. ; Yamaki, T.; Nakanishi, H. *J. Chem. Soc. Chem. Commun.* **1994**, *7*, 855.
- <sup>4</sup> Corey, E.J.; Helal, C. J. Angew. Chem. Int. Ed. 1998, 1986.
- <sup>5</sup> Matos, K.; Corella, J. A.; Burkhardt, E.R.; Nettles, S.M. U.S. 6,218,585.

<sup>6</sup> Amedia Jr., J. C.; Bernard, P.J.; Fountain, M.; Van Wagenen Jr., G. Syn. Commun. **1999**, *29*, 2377.

- <sup>7</sup> Houpis, I.N.; et al. Tetrahedron Letters 1993, 34, 2593. Couturier, M.; et al. Org. Lett. 2001, 3, 465. Couturier, M.; et al. Org. Proc. R&D 2002, 6, 42-46.
- <sup>8</sup> Laaziri, H.; Bromm, L.O.; Lhermitte, F.; Gschwind, R.M.; Knochel, P. *J. Am. Chem. Soc.* **1999**, *121*, 6940.
- <sup>9</sup> Knochel, P. et al. J. Am. Chem. Soc. 2000, 122, 10218.

The information and recommendations contained in this document are presented in good faith and believed to be reliable, but shall not be part of the terms and conditions of sale of any Callery Chemical product. This product may present unknown health hazards and should be used with caution. Although certain hazards are described herein, we cannot guarantee that these are the only hazards that exist. Final determination of the suitability of the product is the sole responsibility of the user. Users of the product should satisfy themselves that the conditions and methods of use assure that the product is used safely. NO REPRESENTATIONS OR WARRANTIES, EITHER EXPRESS OR IMPLIED, OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE OR ANY OTHER NATURE ARE MADE HEREUNDER WITH RESPECT TO THE INFORMATION CONTAINED HEREIN OR THE PRODUCT TO WHICH THE INFORMATION REFERS. Nothing herein is intended as a recommendation to use our products in any manner that would result in the infringement of any patent. We assume no liability for any customer's violation of patent or other rights. The customer should make his own patent investigation relative to his proposed use.